

TO : All Companies

Notice

Subj / Inv No. 64/2022/54

Dear Sirs :

We would like to inform you that technical specifications of a/m equipment have been updated as closing date of the tender has been extended up to 29/11/2022 instead of 15/11/2022 Thanking you for your cooperation .

❖ **Application: middle ear analyzer.**

❖ **FDA clearance or CE MARK(MDD) : Yes FDA Preferred**

CONFIGURATION	Stand alone , or Desk Top
WARM-UP TIME, min	≤ 1
<u>ADMITTANCE MEASUREMENTS</u> Probe tone, Hz	226, 1,000 others to be specified
<u>Accuracy</u> Tympanometry, % Reflex, %	±5 or less ±3 or less
PRESSURE MEASUREMENTS Range, daPa	≥ (+200~-600).
REFLEX MEASUREMENTS Stimuli	Ipsi + Contra (500,1000,2000,4000)Hz
Total harmonic distortion, %	<5
Intensity, dB HL	(35-100) or wider range
TESTS AVAILABLE	<ul style="list-style-type: none">•Automatic and manual tympanometry.•Automatic and manual acoustic reflex (Ipsi/Contra).•Reflex decay (Ipsi/Contra).•Eustachian tube function open &close tympanic membrane.•Acoustic reflex latency.
Quick test	YES
CALIBRATION METHOD	Automatic & Manual
Connectivity	USB OR RS232C
Display	≥ 4" LCD or better
Reports	YES
STANDARD ACCESSORIES	Printer paper audiometry headset (for Contra) ear tips startup Qty ≥ 100 piece Response hand switch on probe.
Input power:	220/240 VAC, 50/60 Hz single phase. Three Pin G type Plug, Voltage stabilizing and over current protection circuit.
Environmental requirements	the equipment suitable for work in the climate conditions in Iraq in terms of temperature & humidity.

B. RQDS

الى : - كافة الشركات

تنوية
م/ المناقصة المرقمه (64/2022/54)
اسم الجهاز Tympanometry

تحية طيبة ...

نود اعلامكم انه تم تحديث المواصفات الفنية للجهاز انفاً و تمديد موعد غلق المناقصة ليكون موعد الغلق ٢٠٢٢/١١/٢٩ بدلا من ٢٠٢٢/١١/١٥ ، شاكرين لكم تعاونكم معنا .

❖ Application: middle ear analyzer.

❖ FDA clearance or CE MARK(MDD) : Yes FDA Preferred

CONFIGURATION	Stand alone , or Desk Top
WARM-UP TIME, min	≤ 1
<u>ADMITTANCE MEASUREMENTS</u> Probe tone, Hz	226, 1,000 others to be specified
<u>Accuracy</u> Tympanometry, % Reflex, %	±5 or less ±3 or less
PRESSURE MEASUREMENTS Range, daPa	≥ (+200~-600).
REFLEX MEASUREMENTS Stimuli	Ipsi + Contra (500,1000,2000,4000)Hz
Total harmonic distortion, %	<5
Intensity, dB HL	(35-100) or wider range
TESTS AVAILABLE	<ul style="list-style-type: none">•Automatic and manual tympanometry.•Automatic and manual acoustic reflex (Ipsi/Contra).•Reflex decay (Ipsi/Contra).•Eustachain tube function open &close tympanic membrane.•Acoustic reflex latency.
Quick test	YES
CALIBRATION METHOD	Automatic & Manual
Connectivity	USB OR RS232C
Display	≥ 4" LCD or better
Reports	YES
STANDARD ACCESSORIES	Printer paper audiometry headset (for Contra) ear tips startup Qty ≥ 100 piece Response hand switch on probe.
Input power:	220/240 VAC, 50/60 Hz single phase. Three Pin G type Plug, Voltage stabilizing and over current protection circuit.
Environmental requirements	the equipment suitable for work in the climate conditions in Iraq in terms of temperature & humidity.

... مع الاحترام